

For Research Use Only

# HURP Recombinant monoclonal antibody, PBS Only

Catalog Number: 86704-1-PBS



## Basic Information

Catalog Number:

86704-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2667

GenBank Accession Number:

BC010658

GeneID (NCBI):

9787

UNIPROT ID:

Q15398

Full Name:

discs, large (Drosophila) homolog-associated protein 5

Calculated MW:

846 aa, 95 kDa

Observed MW:

120 kDa

Purification Method:

Protein A purification

CloneNo.:

251623E3

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86704-1-RR (HURP antibody) at dilution of 1:30000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86704-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86704-1-RR against Human HURP were performed. The affinity constant is 3.84 nM.