

For Research Use Only

GRP75 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86721-3-PBS



Basic Information

Catalog Number: 86721-3-PBS	GenBank Accession Number: BC000478	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 3313	CloneNo.: 251670C11
Source: Rabbit	UNIPROT ID: P38646	
Isotype: IgG	Full Name: heat shock 70kDa protein 9 (mortalin)	
Immunogen Catalog Number: AG6673	Calculated MW: 74 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean GRP75 concentration was determined to be 8.4 ng/mL in Hela cell extract based on a 1.2 mg/mL extract load, 22.6 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02623-1, GRP75 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86721-3-PBS. Detection antibody: 86721-2-PBS. Standard: Ag6673. Range: 0.312-40 ng/mL.