

For Research Use Only

IRG1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 86725-1-PBS



Basic Information

Catalog Number: 86725-1-PBS	GenBank Accession Number: NM_001258406	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 730249	CloneNo.: 251531C12
Isotype: IgG	UNIPROT ID: A6NK06	
Immunogen Catalog Number: AG28913	Full Name: immunoresponsive 1 homolog (mouse)	
	Calculated MW: 53 kDa	
	Observed MW: 53 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human, mouse

Background Information

The host Immune-Responsive Gene 1 (IRG1; also called Acod1) is a mitochondrial enzyme induced under inflammatory conditions that produces the metabolite itaconate by decarboxylating cis-aconitate, a TCA cycle intermediate. Gene expression profiling studies of murine macrophages and microglial cells have revealed that IRG1 is highly expressed under pro-inflammatory conditions. Furthermore, IRG1 is highly expressed in the pregnant uterus during the early events leading to implantation, the specific phase of pregnancy in which high levels of inflammatory cytokines are secreted. IRG1 localizes to the mitochondria and may represent a key link between immunological and metabolic processes. IRG1 has crucial functions in embryonic implantation and neurodegeneration. Also, IRG1 promotes endotoxin tolerance by increasing A20 expression in macrophages via increased ROS production. (PMID: 25640654)

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



LPS treated RAW 264.7 cells were subjected to SDS PAGE followed by western blot with 86725-1-RR (IRG1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86725-1-PBS in a different storage buffer formulation.

Immunofluorescent analysis of (4% PFA) fixed LPS treated RAW 264.7 cells using IRG1 antibody (86725-1-RR, Clone: 251531C12) at dilution of 1:500 and CoraLite@488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 86725-1-PBS in a different storage buffer formulation.