

For Research Use Only

Mouse Cd2ap Recombinant antibody, PBS Only (Detector)

Catalog Number: 86727-2-PBS



Basic Information

Catalog Number:

86727-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0622

GenBank Accession Number:

BC157901

GeneID (NCBI):

12488

UNIPROT ID:

Q9JLQ0

Full Name:

CD2-associated protein

Calculated MW:

637 aa, 70 kDa

Purification Method:

Protein A purification

CloneNo.:

251428G6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Cd2ap concentration was determined to be 948.7 ng/mL in Mouse kidney tissue extract based on a 8.4 mg/mL extract load.



Cytometric bead array standard curve of MP02649-1, MOUSE Cd2ap Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86727-3-PBS. Detection antibody: 86727-2-PBS. Standard: Ag0622. Range: 0.781-100 ng/mL.