

For Research Use Only

TRIM47 Recombinant monoclonal antibody, PBS Only

Catalog Number:86728-1-PBS



Basic Information

Catalog Number: 86728-1-PBS	GenBank Accession Number: BC017299	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 91107	CloneNo.: 251436A2
Isotype: IgG	UNIPROT ID: Q96LD4	
Immunogen Catalog Number: AG22868	Full Name: tripartite motif-containing 47	
	Calculated MW: 638 aa, 70 kDa	
	Observed MW: 70-75 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

TRIM47/GOA is a member of the tripartite motif RBCC/Trim domain family that contains a RING-finger domain, B-boxes, and an alpha-helical coiled-coil region. TRIM47/GOA was identified as a protein overexpressed in astrocytoma tumor cells and is thought to play a role in dedifferentiation associated with tumorigenesis. TRIM47, defined as E3 ubiquitin ligases, physically binds to SMAD4 and enhances the ubiquitination and degradation of SMAD4 (PMID: 39117713, 35513381).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 86728-1-RR (TRIM47 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86728-1-PBS in a different storage buffer formulation.



OVCAR-3 cells were subjected to SDS PAGE followed by western blot with 86728-1-RR (TRIM47 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86728-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 86728-1-RR (TRIM47 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86728-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86728-1-RR against Human TRIM47 were performed. The affinity constant is 1.77 nM.