

For Research Use Only

TRIM47 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86728-2-PBS



Basic Information

| | | |
|---|---|---|
| Catalog Number: 86728-2-PBS | GenBank Accession Number: BC017299 | Purification Method: Protein A purification |
| Source: Rabbit | GeneID (NCBI): 91107 | CloneNo.: 251436B11 |
| Isotype: IgG | UNIPROT ID: Q96LD4 | |
| Immunogen Catalog Number: AG22868 | Full Name: tripartite motif-containing 47 | |
| | Calculated MW: 638 aa, 70 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TRIM47 concentration was determined to be 6.8 ng/mL in U2OS cell extract based on a 1.2 mg/mL extract load, 9.5 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 1.5 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02626-1, TRIM47 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86728-3-PBS. Detection antibody: 86728-2-PBS. Standard: Ag22868. Range: 0.312-40 ng/mL.