

For Research Use Only



Angiotensinogen/AGT Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86729-3-PBS

Basic Information

| | | |
|-------------------------------------|---|--|
| Catalog Number: 86729-3-PBS | GenBank Accession Number: NM_000029.4 | Purification Method: Protein A purification |
| Source: Rabbit | GeneID (NCBI): 183 | CloneNo.: 251419B3 |
| Isotype: IgG | UNIPROT ID: P01019 | |
| Immunogen Catalog Number: EG3099 | Full Name: angiotensinogen (serpin peptidase inhibitor, clade A, member 8) | |
| | Calculated MW: 52 kDa | |

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Serum of six individual healthy human donors was measured. The Angiotensinogen/AGT concentration of detected samples was determined to be 6,000.61 ng/mL with a range of 118.04 - 12,267.72 ng/mL.



Sandwich ELISA standard curve of MP02816-1, Human Angiotensinogen/AGT Recombinant Matched Antibody Pair - PBS only. 86729-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3099. 86729-3-PBS was HRP conjugated as the detection antibody. Range: 0.313-10 ng/mL.