

For Research Use Only

NCF1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:86733-1-PBS



Basic Information

Catalog Number: 86733-1-PBS	GenBank Accession Number: BC002816	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 653361	CloneNo.: 251451A12
Isotype: IgG	UNIPROT ID: P14598	
Immunogen Catalog Number: AG28089	Full Name: neutrophil cytosolic factor 1	
	Calculated MW: 45 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean NCF1/p47 phox concentration was determined to be 133.2 ng/mL in U-937 cell extract based on a 2.6 mg/mL extract load.



Cytometric bead array standard curve of MP02620-1, NCF1/p47 phox Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86733-2-PBS. Detection antibody: 86733-1-PBS. Standard: Ag28089. Range: 1.25-80 ng/mL.