

For Research Use Only

NCF1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86733-2-PBS



Basic Information

Catalog Number:

86733-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28089

GenBank Accession Number:

BC002816

GeneID (NCBI):

653361

UNIPROT ID:

P14598

Full Name:

neutrophil cytosolic factor 1

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

251451C1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean NCF1/p47 phox concentration was determined to be 133.2 ng/mL in U-937 cell extract based on a 2.6 mg/mL extract load.



Cytometric bead array standard curve of MP02620-1, NCF1/p47 phox Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86733-2-PBS. Detection antibody: 86733-1-PBS. Standard: Ag28089. Range: 1.25-80 ng/mL.