

For Research Use Only

NCF1/p47 phox Recombinant monoclonal antibody

Catalog Number:86733-3-RR



Basic Information

Catalog Number:

86733-3-RR

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28089

GenBank Accession Number:

BC002816

GeneID (NCBI):

653361

UNIPROT ID:

P14598

Full Name:

neutrophil cytosolic factor 1

Calculated MW:

45 kDa

Observed MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

252310G3

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : Ramos cells, Daudi cells, Raji cells, RAW 264.7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86733-3-RR (NCF1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.