

For Research Use Only

# Di-Methyl-Histone H3 (Lys79) Recombinant monoclonal antibody

Catalog Number:86734-1-RR



## Basic Information

Catalog Number:

86734-1-RR

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC066245

GeneID (NCBI):

8350

UNIPROT ID:

P68431

Full Name:

histone cluster 1, H3a

Observed MW:

15 kDa

Purification Method:

Protein A purification

CloneNo.:

251570E11

Recommended Dilutions:

WB: 1:500-1:2000

Dot Blot: 1:10-1:100

ChIP-qPCR: 1:10-1:100

## Applications

Tested Applications:

WB, Dot Blot, ELISA, ChIP-qPCR

Species Specificity:

human, mouse

Positive Controls:

WB : mouse liver tissue, MCF-7 cells

Dot Blot: //,

ChIP-qPCR : HeLa cells,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data




Dot blot analysis was used to confirm the specificity of 86734-1-RR Di-Methyl-Histone H3 (Lys79) antibody. peptides were spotted onto NC and probed with antibody at 1 µg/ml. The amount of peptide ( µ g/mL) spotted is indicated next to each row.



Various lysates were subjected to SDS PAGE followed by western blot with 86734-1-RR (Di-Methyl-Histone H3 (Lys79) antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Chromatin was prepared from HeLa cells. Cells were fixed with formaldehyde for 10 minutes. The CHIP was performed with 15 µg of cross-linked chromatin, 5 µg of Di-Methyl-Histone H3 (Lys79) (86734-1-RR) or 5 µg of Normal Rabbit IgG (98136-1-RR), and 20 µl of Protein A Magarose Beads. The immunoprecipitated DNA was quantified by real-time PCR.



mouse liver tissue were subjected to SDS PAGE followed by western blot with 86734-1-RR (Di-Methyl-Histone H3 (Lys79) antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.