

For Research Use Only

# PHPT1 Recombinant monoclonal antibody, PBS Only

Catalog Number:86736-3-PBS



## Basic Information

Catalog Number:

86736-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3778

GenBank Accession Number:

BC024648

GeneID (NCBI):

29085

UNIPROT ID:

Q9NRX4

Full Name:

phosphohistidine phosphatase 1

Calculated MW:

125 aa, 14 kDa

Observed MW:

14 kDa

Purification Method:

Protein A purification

CloneNo.:

251432C1

## Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86736-3-RR (PHPT1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86736-3-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HEK-293 cells using PHPT1 antibody (86736-3-RR, Clone: 251432C1) at dilution of 1:200 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 86736-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86736-3-RR against Human PHPT1 were performed. The affinity constant is 2.39 nM.