

For Research Use Only

RNF43 Recombinant monoclonal antibody, PBS Only

Catalog Number: 86740-1-PBS



Basic Information

Catalog Number:

86740-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG6490

GenBank Accession Number:

NM_001305544.1

GeneID (NCBI):

54894

UNIPROT ID:

Q68DV7

Full Name:

ring finger protein 43

Calculated MW:

86 kDa

Observed MW:

86 kDa

Purification Method:

Protein A purification

CloneNo.:

251708C1

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86740-1-RR (RNF43 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86740-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HEK-293 cells using RNF43 antibody (86740-1-RR, Clone: 251708C1) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 86740-1-PBS in a different storage buffer formulation.