

For Research Use Only

PDPN Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86750-5-PBS



Basic Information

Catalog Number:

86750-5-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5024

GenBank Accession Number:

BC014668

GeneID (NCBI):

10630

UNIPROT ID:

Q86YL7-1

Full Name:

podoplanin

Purification Method:

Protein A purification

CloneNo.:

251406C8

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Plasma of eight individual healthy human donors was measured. The Podoplanin concentration of detected samples was determined to be 85.9 pg/mL with a range of 37.1-122.1 pg/mL.



Sandwich ELISA standard curve of MP02738-1, Human Podoplanin Recombinant Matched Antibody Pair - PBS only. 86750-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg5024. 86750-5-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL.