

For Research Use Only

Serpin A4/Kallistatin Recombinant monoclonal antibody

Catalog Number:86758-1-RR



Basic Information

Catalog Number: 86758-1-RR	GenBank Accession Number: NM_001289033.1	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 5267	CloneNo.: 251716B4
Isotype: IgG	UNIPROT ID: P29622	Recommended Dilutions: WB: 1:5000-1:50000
Immunogen Catalog Number: EG4707	Full Name: serpin peptidase inhibitor, clade A (alpha-1 antitrypsin, antitrypsin), member 4	
	Calculated MW: 49 kDa	
	Observed MW: 49 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human plasma, mouse plasma
Species Specificity: human, mouse	

Background Information

Serpin A4 is also named Kallistatin and Peptidase inhibitor 4 (PI-4). Kallistatin was first identified in human plasma as a tissue kallikrein-binding protein and a serine proteinase inhibitor. Kallistatin via its two structural elements regulates differential signaling cascades, and thus a wide spectrum of biological functions (PMID: 28427524). Kallistatin via the active site promotes neovascularization by stimulating VEGF levels in endothelial progenitor cells. Second, kallistatin inhibits or induces apoptosis depending on cell types. Kallistatin via its heparin-binding site inhibits angiogenesis by blocking vascular endothelial growth factor (VEGF)-induced growth, migration and adhesion of endothelial cells. Kallistatin can act as an antioxidant or pro-oxidant (PMID: 28742513).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86758-1-RR (Serpin A4/Kallistatin antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.