

For Research Use Only

Mouse FABP1 Recombinant antibody, PBS Only (Capture)

Catalog Number:86759-2-PBS



Basic Information

Catalog Number: 86759-2-PBS	GenBank Accession Number: NM_017399.4	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 14080	CloneNo.: 251790E6
Source: Rabbit	UNIPROT ID: P12710	
Isotype: IgG	Full Name: fatty acid binding protein 1, liver	
Immunogen Catalog Number: EG7278	Calculated MW: 14 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FABP1 concentration was determined to be 8.4 ng/mL in mouse serum.



Cytometric bead array standard curve of MP02652-1, MOUSE FABP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86759-2-PBS. Detection antibody: 86759-1-PBS. Standard: Eg7278. Range: 0.312-40 ng/mL.