

For Research Use Only

FABP1 Recombinant monoclonal antibody

Catalog Number:86759-3-RR



Basic Information

Catalog Number: 86759-3-RR	GenBank Accession Number: NM_017399.4	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 14080	CloneNo.: 251790G4
Isotype: IgG	UNIPROT ID: P12710	Recommended Dilutions: WB: 1:5000-1:50000
Immunogen Catalog Number: EG7278	Full Name: fatty acid binding protein 1, liver	
	Calculated MW: 14 kDa	
	Observed MW: 14 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse liver tissue, rat liver tissue, pig liver tissue, human kidney tissue
Species Specificity: human, mouse, rat, pig	

Background Information

FABP1, liver fatty acid-binding protein, is abundant in cytoplasm that regulates lipid transport and metabolism. It plays a role in lipoprotein-mediated cholesterol uptake in hepatocytes (PMID:25732850). FABP1 is mainly expressed in hepatocytes, enterocytes and to a lesser degree in renal tubular cells, associated with liver injury (PMID: 15653098). Levels of FABP1 have been found to be elevated in patients with hepatocyte injury secondary to alcohol or drug toxicity (PMID: 14563446).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86759-3-RR (FABP1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.