

For Research Use Only

Mouse RANK/TNFRSF11A Recombinant monoclonal antibody, PBS Only (Detector)



Catalog Number: 86769-1-PBS

Basic Information

Catalog Number: 86769-1-PBS	GenBank Accession Number: NM_009399.3	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 21934	CloneNo.: 250104A4
Source: Rabbit	UNIPROT ID: O35305	
Isotype: IgG	Full Name: tumor necrosis factor receptor superfamily, member 11a	
Immunogen Catalog Number: EG1829	Calculated MW: 67kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean RANK/TNFRSF11A concentration was determined to be 5.88 ng/mL in RAW 264.7 cell extract based on a 1.30 mg/mL extract load.



Sandwich ELISA standard curve of MP02637-1, Mouse RANK/TNFRSF11A Recombinant Matched Antibody Pair - PBS only. 86769-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg1829. 86769-1-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL.