

For Research Use Only

TBR1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:86772-2-PBS



Basic Information

Catalog Number: 86772-2-PBS	GenBank Accession Number: BC104844	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 10716	CloneNo.: 25143484
Isotype: IgG	UNIPROT ID: Q16650	
Immunogen Catalog Number: AG14935	Full Name: T-box, brain, 1	
	Calculated MW: 682 aa, 74 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TBR1 concentration was determined to be 2.9 ng/mL in SH-SY5Y cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02778-1, Human TBR1 Recombinant Matched Antibody Pair - PBS only. 86772-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag14935. 86772-2-PBS was HRP conjugated as the detection antibody. Range: 31.3-2000 pg/mL.