

For Research Use Only

HFE2/Hemojuvelin Recombinant monoclonal antibody, PBS Only

Catalog Number:86775-3-PBS



Basic Information

Catalog Number: 86775-3-PBS	GenBank Accession Number: BC017926	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 148738	CloneNo.: 251477G9
Isotype: IgG	UNIPROT ID: Q6ZVN8	
Immunogen Catalog Number: AG2354	Full Name: hemochromatosis type 2 (juvenile)	
	Calculated MW: 426 aa, 45 kDa	
	Observed MW: 45 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human, mouse

Background Information

HFE2, also known as Hemojuvelin (HJV), is a GPI-anchored protein that binds to new proteins and bone morphogenetic proteins (BMP2 and BMP4) and can be released from the cell membrane. HFE2 plays a key role in the regulation of hepcidin, which regulates the absorption and recycling of iron.(PMID: 17331953; PMID: 17938254)

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86775-3-RR (HFE2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86775-3-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using HFE2 antibody (86775-3-RR, Clone: 251477G9) at dilution of 1:200 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 86775-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86775-3-RR against Human HFE2 were performed. The affinity constant is below 1 pM.