

For Research Use Only

ACTRT2 Recombinant monoclonal antibody, PBS Only

Catalog Number:86782-1-PBS



Basic Information

Catalog Number:

86782-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10597

GenBank Accession Number:

BC029499

GeneID (NCBI):

140625

UNIPROT ID:

Q8TDY3

Full Name:

actin-related protein T2

Calculated MW:

376 aa, 42 kDa

Observed MW:

40-42 kDa

Purification Method:

Protein A purification

CloneNo.:

251447B12

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86782-1-RR (ACTRT2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86782-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86782-1-RR against Human ACTRT2 were performed. The affinity constant is below 1 pM.