

For Research Use Only

PHLDA2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:86784-1-PBS



Basic Information

Catalog Number: 86784-1-PBS	GenBank Accession Number: BC005034	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 7262	CloneNo.: 251671G5
Isotype: IgG	UNIPROT ID: Q53GA4	
Immunogen Catalog Number: AG6331	Full Name: pleckstrin homology-like domain, family A, member 2	
	Calculated MW: 17 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PHLDA2 concentration was determined to be 2.7 ng/mL in L02 cell extract based on a 1.3 mg/mL extract load, 1.4 ng/mL in PC-3 cell extract based on a 1.3 mg/mL extract load.



Cytometric bead array standard curve of MP02666-1, PHLDA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86784-2-PBS. Detection antibody: 86784-1-PBS. Standard: Ag6331. Range: 0.156-20 ng/mL.