

For Research Use Only

# PARD6B Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86786-3-PBS



## Basic Information

Catalog Number:

86786-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5083

GenBank Accession Number:

BC060847

GeneID (NCBI):

84612

UNIPROT ID:

Q9BYG5

Full Name:

par-6 partitioning defective 6 homolog beta (C. elegans)

Calculated MW:

41 kDa

Purification Method:

Protein A purification

CloneNo.:

251629A5

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean PARD6B concentration was determined to be 60.12 ng/mL in HepG2 cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP02746-1, Human PARD6B Recombinant Matched Antibody Pair - PBS only. 86786-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag5083. 86786-3-PBS was HRP conjugated as the detection antibody. Range: 1.56-100 ng/mL