

For Research Use Only

PARD6B Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86786-4-PBS



Basic Information

Catalog Number:

86786-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5083

GenBank Accession Number:

BC060847

GeneID (NCBI):

84612

UNIPROT ID:

Q9BYG5

Full Name:

par-6 partitioning defective 6 homolog beta (C. elegans)

Calculated MW:

41 kDa

Purification Method:

Protein A purification

CloneNo.:

251629A11

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PARD6B concentration was determined to be 60.12 ng/mL in HepG2 cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP02746-1, Human PARD6B Recombinant Matched Antibody Pair - PBS only. 86786-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag5083. 86786-3-PBS was HRP conjugated as the detection antibody. Range: 1.56-100 ng/mL