

For Research Use Only

EPHX2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86790-2-PBS



Basic Information

Catalog Number:

86790-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1283

GenBank Accession Number:

BC013874

GeneID (NCBI):

2053

UNIPROT ID:

P34913

Full Name:

epoxide hydrolase 2, cytoplasmic

Calculated MW:

63 kDa

Purification Method:

Protein A purification

CloneNo.:

251753F5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean EPHX2 concentration was determined to be 15.72 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 24.32 ng/mL in HEK-293T cell extract based on a 1.3 mg/mL extract load and 11.98 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02802-1, Human EPHX2 Recombinant Matched Antibody Pair - PBS only. 86790-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag1283. 86790-2-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL