

For Research Use Only

# FABP1 Recombinant monoclonal antibody

Catalog Number:86792-1-RR



## Basic Information

<b>Catalog Number:</b> 86792-1-RR	<b>GenBank Accession Number:</b> NM_001443.3	<b>Purification Method:</b> Protein A purification
<b>Source:</b> Rabbit	<b>GeneID (NCBI):</b> 2168	<b>CloneNo.:</b> 251782A3
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> P07148	<b>Recommended Dilutions:</b> WB: 1:5000-1:50000
<b>Immunogen Catalog Number:</b> EG7273	<b>Full Name:</b> fatty acid binding protein 1, liver	
	<b>Calculated MW:</b> 14 kDa	
	<b>Observed MW:</b> 14 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : HepG2 cells, mouse liver tissue
<b>Species Specificity:</b> human, mouse	

## Background Information

FABP1, liver fatty acid-binding protein, is abundant in cytoplasm that regulates lipid transport and metabolism. It plays a role in lipoprotein-mediated cholesterol uptake in hepatocytes (PMID:25732850). FABP1 is mainly expressed in hepatocytes, enterocytes and to a lesser degree in renal tubular cells, associated with liver injury (PMID: 15653098). Levels of FABP1 have been found to be elevated in patients with hepatocyte injury secondary to alcohol or drug toxicity (PMID: 14563446).

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol, pH7.3  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86792-1-RR (FABP1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.