

For Research Use Only

AP3S2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:86795-2-PBS



Basic Information

Catalog Number:

86795-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7565

GenBank Accession Number:

BC002785

GeneID (NCBI):

10239

UNIPROT ID:

P59780

Full Name:

adaptor-related protein complex 3,
sigma 2 subunit

Calculated MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

251741A7

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean AP3S2 concentration was determined to be 9.81 ng/mL in mouse brain tissue extract based on a 2.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02806-1, Human AP3S2 Recombinant Matched Antibody Pair - PBS only. 86795-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag7565. 86795-1-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL