

For Research Use Only

AZI1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:86809-1-PBS



Basic Information

Catalog Number: 86809-1-PBS	GenBank Accession Number: BC136578	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 22994	CloneNo.: 251839B2
Source: Rabbit	UNIPROT ID: Q9UPN4	
Isotype: IgG	Full Name: 5-azacytidine induced 1	
Immunogen Catalog Number: AG22699	Calculated MW: 1083 aa, 122 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean AZI1 concentration was determined to be 13.4 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02687-1, AZI1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86809-2-PBS. Detection antibody: 86809-1-PBS. Standard: Ag22699. Range: 0.625-80 ng/mL.