

For Research Use Only

LTBR Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86823-2-PBS



Basic Information

Catalog Number: 86823-2-PBS	GenBank Accession Number: BC026262	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 4055	CloneNo.: 251739F7
Source: Rabbit	UNIPROT ID: P36941	
Isotype: IgG	Full Name: lymphotoxin beta receptor (TNFR superfamily, member 3)	
Immunogen Catalog Number: EG5126	Calculated MW: 435 aa, 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean LTBR concentration was determined to be 72.3 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 35.0 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02693-1, LTBR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86823-3-PBS. Detection antibody: 86823-2-PBS. Standard: Eg5126. Range: 0.781-100 ng/mL