

For Research Use Only

LTBR Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86823-4-PBS



Basic Information

Catalog Number:

86823-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5126

GenBank Accession Number:

BC026262

GeneID (NCBI):

4055

UNIPROT ID:

P36941

Full Name:

lymphotoxin beta receptor (TNFR superfamily, member 3)

Calculated MW:

435 aa, 47 kDa

Purification Method:

Protein A purification

CloneNo.:

251739F4

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean LTBR concentration was determined to be 14.22 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load and 17.35 ng/mL in Caki-1 cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02693-2, Human LTBR Recombinant Matched Antibody Pair - PBS only. 86823-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg5126. 86823-4-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL