

For Research Use Only

TXNDC9 Recombinant monoclonal antibody, PBS Only

Catalog Number:86826-5-PBS



Basic Information

Catalog Number: 86826-5-PBS	GenBank Accession Number: BC005968	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 10190	CloneNo.: 251866B11
Isotype: IgG	UNIPROT ID: O14530	
Immunogen Catalog Number: AG19509	Full Name: thioredoxin domain containing 9	
	Calculated MW: 226 aa, 27 kDa	
	Observed MW: 25-27 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

TXNDC9, also known as APACD, has a negative effect on protein folding by significantly reducing the activity of chaperone protein TCP1 complex ATPase activity. Including actin or tubulin.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86826-5-RR (TXNDC9 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86826-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86826-5-RR against Human TXNDC9 were performed. The affinity constant is below 1 pM.