

For Research Use Only

RAB34 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:86839-2-PBS



Basic Information

Catalog Number:

86839-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26696

GenBank Accession Number:

BC016841

GeneID (NCBI):

83871

UNIPROT ID:

Q9BZG1

Full Name:

RAB34, member RAS oncogene family

Calculated MW:

251 aa, 28 kDa

Purification Method:

Protein A purification

CloneNo.:

251867F11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean RAB34 concentration was determined to be 30.5 ng/mL in A431 cell extract based on a 1.3 mg/mL extract load, 18.7 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 39.6 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02698-1, RAB34 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86839-2-PBS. Detection antibody: 86839-1-PBS. Standard: Ag26696. Range: 0.469-60 ng/mL.