

For Research Use Only

BUB1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86843-5-PBS



Basic Information

Catalog Number:

86843-5-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4138

GenBank Accession Number:

BC028201

GeneID (NCBI):

699

UNIPROT ID:

O43683

Full Name:

budding uninhibited by
benzimidazoles 1 homolog (yeast)

Calculated MW:

1085 aa, 122 kDa

Purification Method:

Protein A purification

CloneNo.:

251580D11

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean BUB1 concentration was determined to be 304.97 ng/mL in MKN-45 cell extract based on a 2.0 mg/mL extract load.



Sandwich ELISA standard curve of MP02812-1, Human BUB1 Recombinant Matched Antibody Pair - PBS only. 86843-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag4138. 86843-4-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL