

For Research Use Only

# SPIN1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86848-2-PBS



## Basic Information

Catalog Number: 86848-2-PBS	GenBank Accession Number: BC013571	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 10927	CloneNo.: 251836G10
Isotype: IgG	UNIPROT ID: Q9Y657	
Immunogen Catalog Number: AG2750	Full Name: spindlin 1	
	Calculated MW: 262 aa, 30 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean SPIN1 concentration was determined to be 23.1 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02730-1, SPIN1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86848-3-PBS. Detection antibody: 86848-2-PBS. Standard: Ag2750. Range: 0.469-60 ng/mL.