

For Research Use Only

HBG1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 86851-1-PBS



Basic Information

Catalog Number:

86851-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22776

GenBank Accession Number:

BC010913

GeneID (NCBI):

3047

UNIPROT ID:

P69891

Full Name:

hemoglobin, gamma A

Calculated MW:

147 aa, 16 kDa

Observed MW:

14-16 kDa

Purification Method:

Protein A purification

CloneNo.:

251754H5

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

HBG1, also named as Hb F Agamma, HBGA, HBGR and HSGGL, belongs to the globin family. HBG1 makes up the fetal hemoglobin F, in combination with alpha chains. Some gamma variants can cause severe jaundice and cyanosis in premature and new born babies. The antibody is specific to HBG1 and HBG2. The antibody has no cross reaction to HBB and HBD. HBG1 protein migrates around 16 kDa. Heterodimer of 32 kDa may also be observed. This antibody could identify HBG1 and HBG2.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86851-1-RR (HBG1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86851-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed K-562 cells using HBG1 antibody (86851-1-RR, Clone: 251754H5) at dilution of 1:200 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 86851-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86851-1-RR against Human HBG1 were performed. The affinity constant is 0.941 nM.