

For Research Use Only

FBP17 Recombinant monoclonal antibody, PBS Only

Catalog Number:86852-2-PBS



Basic Information

Catalog Number: 86852-2-PBS	GenBank Accession Number: BC101755	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 23048	CloneNo.: 251876G12
Isotype: IgG	UNIPROT ID: Q96RU3	
Immunogen Catalog Number: AG25784	Full Name: formin binding protein 1	
	Observed MW: 80-85 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

FBP17 (formin-binding protein 17) is a new dynamin-associating molecule consisting of several functional domains, including an N-terminal EFC domain and an SH3 domain at the C-terminus. FBP17 is involved in dynamin-mediated endocytosis and has been reported to be required for invadopodia formation and tumor cell invasion.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86852-2-RR (FNBP1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86852-2-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed Jurkat cells using FBP17 antibody (86852-2-RR, Clone: 251876G12) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 86852-2-PBS in a different storage buffer formulation.