

For Research Use Only

TRAIL R3/TNFRSF10C Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86886-2-PBS



Basic Information

Catalog Number: 86886-2-PBS	GenBank Accession Number: NM_003841.5	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 8794	CloneNo.: 251096B3
Isotype: IgG	UNIPROT ID: O14798	
Immunogen Catalog Number: EG4677	Full Name: tumor necrosis factor receptor superfamily, member 10c, decoy without an intracellular domain	
	Calculated MW: 27 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean DcR1/TNFRSF10C concentration was determined to be 2.6 ng/mL in JAR cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02783-1, Human DcR1/TNFRSF10C Recombinant Matched Antibody Pair - PBS only. 86886-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg4677. 86886-2-PBS was HRP conjugated as the detection antibody. Range: 15.6-1000 pg/mL.