

For Research Use Only

FZD7 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86904-2-PBS



Basic Information

Catalog Number:

86904-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9446

GenBank Accession Number:

BC015915

GeneID (NCBI):

8324

UNIPROT ID:

O75084

Full Name:

frizzled homolog 7 (Drosophila)

Calculated MW:

574 aa, 64 kDa

Purification Method:

Protein A purification

CloneNo.:

250969B10

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Frizzled 7 concentration was determined to be 15.12 ng/mL in MDA-MB-231 cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP02779-1, Human Frizzled 7 Recombinant Matched Antibody Pair - PBS only. 86904-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag9446. 86904-2-PBS was HRP conjugated as the detection antibody. Range: 0.469-30 ng/mL.