

For Research Use Only

CGREF1 Recombinant monoclonal antibody, PBS Only

Catalog Number:86905-1-PBS



Basic Information

Catalog Number: 86905-1-PBS	GenBank Accession Number: BC034764	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 10669	CloneNo.: 251879C5
Isotype: IgG	UNIPROT ID: Q99674	
Immunogen Catalog Number: AG35664	Full Name: cell growth regulator with EF-hand domain 1	
	Calculated MW: 301 aa, 32 kDa	
	Observed MW: 20-22 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

CGREF1, also known as CGR11, is predicted to enable calcium ion binding activity and to be involved in negative regulation of cell population proliferation. Human GREF1 protein has 4 clustered 17-amino acid repeats within the C terminus, suggesting a structure resembling 4 alpha-helical domains interrupted by distinct turns. Northern blot analysis of rat tissues demonstrated expression of CGR11 predominantly in the whole brain and kidney, with limited expression in heart, lung, liver, and skeletal muscle and no expression in spleen and testis.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86905-1-RR (CGREF1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86905-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86905-1-RR against Human CGREF1 were performed. The affinity constant is 4.15 nM.