

For Research Use Only

Fibrinogen alpha chain Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86913-1-PBS



Basic Information

Catalog Number: 86913-1-PBS	GenBank Accession Number: NM_000508.4	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 2243	CloneNo.: 251931C8
Isotype: IgG	UNIPROT ID: P02671-1	
Immunogen Catalog Number: EG4643	Full Name: fibrinogen alpha chain	
	Calculated MW: 95kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean Fibrinogen Alpha Chain concentration was determined to be 8.2 ng/mL in Human placenta tissue extract based on a 3.2 mg/mL extract load.



The mean Fibrinogen Alpha Chain concentration was determined to be 455.5 ng/mL in human serum.



Cytometric bead array standard curve of MP02768-1, Fibrinogen Alpha Chain Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86913-2-PBS. Detection antibody: 86913-1-PBS. Standard: Eg4643. Range: 0.625-40 ng/mL.