

For Research Use Only

# KAT2B/PCAF Recombinant monoclonal antibody, PBS Only

Catalog Number: 86917-1-PBS



## Basic Information

<b>Catalog Number:</b> 86917-1-PBS	<b>GenBank Accession Number:</b> BC060823	<b>Purification Method:</b> Protein A purification
<b>Source:</b> Rabbit	<b>GeneID (NCBI):</b> 8850	<b>CloneNo.:</b> 251940G10
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> Q92831	
<b>Immunogen Catalog Number:</b> AG30238	<b>Full Name:</b> K(lysine) acetyltransferase 2B	
	<b>Calculated MW:</b> 93 kDa	
	<b>Observed MW:</b> 93 kDa	

## Applications

**Tested Applications:**  
WB, Indirect ELISA

**Species Specificity:**  
human

## Background Information

KAT2B, also known as PCAF (p300/CBP-associated factor), is a crucial transcriptional coactivator that possesses intrinsic histone acetyltransferase (HAT) activity. It primarily acetylates histone H3 at lysine 9 (H3K9ac) and lysine 14 (H3K14ac), a modification associated with open chromatin and active gene transcription. Functionally, PCAF is recruited to promoters by interacting with various transcription factors (such as p53, NF- $\kappa$ B, and nuclear receptors) and larger complexes like the p300/CBP cointegrator complex. It plays pivotal roles in diverse cellular processes, including cell cycle progression, differentiation, metabolism, DNA damage response, and tumor suppression or promotion in a context-dependent manner. Due to its involvement in critical pathways, PCAF is a significant protein of interest in epigenetics, cancer biology, and neurodevelopmental research.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86917-1-RR (KAT2B antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86917-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86917-1-RR against Human KAT2B/PCAF were performed. The affinity constant is 2.06 nM.