

For Research Use Only

HNRNPUL1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86926-1-PBS



Basic Information

Catalog Number:

86926-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0899

GenBank Accession Number:

BC002564

GeneID (NCBI):

11100

UNIPROT ID:

Q9BUJ2

Full Name:

heterogeneous nuclear
ribonucleoprotein U-like 1

Calculated MW:

120 kDa

Purification Method:

Protein A purification

CloneNo.:

252010C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean HNRNPUL1 concentration was determined to be 87.9 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, 95.9 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 33.9 ng/mL in HuH-7 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02774-1, HNRNPUL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86926-2-PBS. Detection antibody: 86926-1-PBS. Standard: Ag0899. Range: 0.312-40 ng/mL.