

For Research Use Only

IL-16 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86927-1-PBS



Basic Information

Catalog Number:

86927-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4831

GenBank Accession Number:

AAC12732.1

GeneID (NCBI):

3603

Full Name:

interleukin 16 (lymphocyte chemoattractant factor)

Calculated MW:

9 kDa

Purification Method:

Protein A purification

CloneNo.:

251661B4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Human peripheral blood mononuclear cells (PBMC) were cultured unstimulated or stimulated with 10 μ g/mL PHA for 5 days. The mean IL-16 concentration was determined to be 5.5 ng/mL in unstimulated PBMC supernatant, 15.5 ng/mL in PHA stimulated PBMC supernatant.



Cytometric bead array standard curve of MP02772-1, IL-16 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86927-2-PBS. Detection antibody: 86927-1-PBS. Standard: Eg4831. Range: 0.781-100 ng/mL