

For Research Use Only

IL-16 Recombinant monoclonal antibody

Catalog Number: 86927-2-RR



Basic Information

Catalog Number:

86927-2-RR

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4831

GenBank Accession Number:

AAC12732.1

GeneID (NCBI):

3603

Full Name:

interleukin 16 (lymphocyte chemoattractant factor)

Calculated MW:

9 kDa

Observed MW:

80 kDa

Purification Method:

Protein A purification

CloneNo.:

251661G11

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Recombinant protein, Jurkat cells, HuT 78 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86927-2-RR (IL-16 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 86927-2-RR (IL-16 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.