

For Research Use Only

RCOR1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 86954-1-PBS



Basic Information

Catalog Number: 86954-1-PBS	GenBank Accession Number: NM_015156	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 23186	CloneNo.: 252115G9
Isotype: IgG	UNIPROT ID: Q9UKLO	
Immunogen Catalog Number: AG26770	Full Name: REST corepressor 1	
	Calculated MW: 53 kDa	
	Observed MW: 60 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

RCOR1, also named as REST corepressor 1, is a 485 amino acid protein, which belongs to the CoREST family. RCOR1 is an essential component of the BHC complex, a corepressor complex that represses transcription of neuron-specific genes in non-neuronal cells. The BHC complex is recruited at RE1/NRSE sites by REST and acts by deacetylating and demethylating specific sites on histones, thereby acting as a chromatin modifier. The calculated molecular weight of RCOR1 is 53 kDa, but the phosphorylated RCOR1 is about 60-65 kDa.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com


This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86954-1-RR (RCOR1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86954-1-PBS in a different storage buffer formulation.

Immunofluorescent analysis of (4% PFA) fixed HeLa cells using RCOR1 antibody (86954-1-RR, Clone: 252115G9) at dilution of 1:1000 and Coralite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 86954-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86954-1-RR against Human RCOR1 were performed. The affinity constant is 1.20 nM.