

For Research Use Only

Filamin A Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86959-2-PBS



Basic Information

Catalog Number: 86959-2-PBS	GenBank Accession Number: BC014654	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 2316	CloneNo.: 252084C6
Isotype: IgG	UNIPROT ID: P21333	
Immunogen Catalog Number: AG1581	Full Name: filamin A, alpha (actin binding protein 280)	
	Calculated MW: 2647 aa, 280 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean Filamin A concentration was determined to be 321.5 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.



The mean Filamin A concentration was determined to be 115.6 ng/mL in human plasma.



Cytometric bead array standard curve of MP02793-1, Filamin A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86959-2-PBS. Detection antibody: 86959-1-PBS. Standard: Ag1581. Range: 0.781-100 ng/mL