

For Research Use Only

Filamin A Recombinant monoclonal antibody, PBS Only

Catalog Number:86959-3-PBS



Basic Information

Catalog Number: 86959-3-PBS	GenBank Accession Number: BC014654	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 2316	CloneNo.: 252084C12
Isotype: IgG	UNIPROT ID: P21333	
Immunogen Catalog Number: AG1581	Full Name: filamin A, alpha (actin binding protein 280)	
	Calculated MW: 2647 aa, 280 kDa	
	Observed MW: 280 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human, mouse

Background Information

FLNA encodes a 280 kDa widely expressed actin binding protein called filamin A which can crosslinks actin filaments and links actin filaments to membrane glycoproteins to form a three-dimensional network (PMID:29931263). And filamin A could interact with many other proteins, involved in receptor activation, cell migration and adhesion, cell proliferation, inflammation and tumorigenesis, studies have been reported that FLNA overexpressed in multiple types of tumors, suggesting FLNA may involved in cancer aggressiveness(PMID:29100390). What's more, FLNA is a causative gene of periventricular nodular heterotopia (PNH) (PMID: 29062687)

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86959-3-RR (FLNA antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86959-3-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed HeLa cells using Filamin A antibody (86959-3-RR, Clone: 252084C12) at dilution of 1:1000 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 86959-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86959-3-RR against Human Filamin A were performed. The affinity constant is 6.35 nM.