

For Research Use Only

KAT6B/MORF Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:86962-8-PBS



Basic Information

Catalog Number:

86962-8-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29842

GenBank Accession Number:

BC150618

GeneID (NCBI):

23522

UNIPROT ID:

Q8WYB5

Full Name:

MYST histone acetyltransferase
(monocytic leukemia) 4

Calculated MW:

2073 aa, 231 kDa

Purification Method:

Protein A purification

CloneNo.:

251928G7

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean KAT6B/MORF concentration was determined to be 261.93 ng/mL in MCF-7 cell extract based on a 1.20 mg/mL extract load and 224.0 ng/mL in MDA-MB-231 cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP02869-1, Human KAT6B/MORF Recombinant Matched Antibody Pair - PBS only. 86962-8-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag29842. 86962-7-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL.