

For Research Use Only

# MYST3 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86963-5-PBS



## Basic Information

Catalog Number:

86963-5-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29729

GenBank Accession Number:

BC142659

GeneID (NCBI):

7994

UNIPROT ID:

Q92794

Full Name:

MYST histone acetyltransferase  
(monocytic leukemia) 3

Calculated MW:

2004 aa, 225 kDa

Purification Method:

Protein A purification

CloneNo.:

251927H8

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean MYST3 concentration was determined to be 3.35 ng/mL in HT-1080 cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP02781-1, Human MYST3 Recombinant Matched Antibody Pair - PBS only. 86963-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag29729. 86963-5-PBS was HRP conjugated as the detection antibody. Range: 0.313-10 ng/mL