

For Research Use Only

CNTN3 Recombinant monoclonal antibody, PBS Only

Catalog Number:86965-1-PBS



Basic Information

Catalog Number:

86965-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5928

GenBank Accession Number:

NM_020872.2

GeneID (NCBI):

5067

UNIPROT ID:

Q9P232

Full Name:

contactin 3 (plasmacytoma associated)

Calculated MW:

113 kDa

Observed MW:

140 kDa

Purification Method:

Protein A purification

CloneNo.:

251948A2

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat, pig

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86965-1-RR (CNTN3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86965-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 86965-1-RR (CNTN3 antibody) at dilution of 1:500 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 86965-1-PBS in a different storage buffer formulation.